Scan for more information

Advazorb®/Advazorb® Lite

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Advanc**Medical**

Advazorb®/Advazorb® Lite

Hydrophilic foam dressing with film backing

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
Product description

A soft, conforming non-adherent hydrophilic polyurethane foam wound dressing with pink polyurethane film backing. Advazorb® has high fluid absorbency and retention characteristics.

Indications

Suitable for a wide range of exuding wounds, granulating wounds, leg ulcers, pressure sores and for protection of bony prominences.

Contra-indications

Arterial bleeds and heavily bleeding wounds.

Directions for use

Place directly to the wound surface pink film side up, ensuring the absorbent pad covers the entire wound area with an overlap of 1-2cm. Secure in place with tape, appropriate bandage or film dressing. Can be used under compression.

Frequency of change

Can be left in place for up to 7 days but should be changed immediately if strike through is apparent or when dressing reaches capacity. If exudate is visible around the edges of the dressing this is a clear indication that a dressing change is required. Clinical observation is necessary to determine required frequency of change.

Warnings and Precautions

Advazorb® is a single use product. If used on more than one patient, cross-contamination or infection may result. Opening the dressing pack compromises the sterile barrier. Any unused dressings should be discarded.

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