

Scan for more information



Advazorb[®] Heel

Absorbent foam dressing with film backing

Instructions for use



Advancis Medical, Lowmoor Business Park,
Kirkby-in-Ashfield, Nottingham, NG17 7JZ, England

Tel: +44 (0) 1623 751500

Fax: +44 (0) 871 264 8238

Email: info@advancis.co.uk

Web: www.advancis.co.uk



Advazorb® Heel

Absorbent foam dressing with film backing

EN



Product Description

Advazorb® Heel is a conformable, non-adherent, absorbent, polyurethane foam dressing with a breathable film backing designed to fit around the heel. A non-adherent foam dressing with excellent fluid handling properties. Ideal for use under compression bandaging with its fluid retention, low profile and low friction film backing. Advazorb® Heel also contributes to a pressure relieving protocol when used in conjunction with pressure relieving devices.

Indications

As a wound dressing, Advazorb® Heel is suitable for acute and chronic wounds including: cuts and abrasions, superficial burns, surgical wounds, leg ulcers, pressure ulcers and diabetic ulcers. Advazorb® Heel can be used as part of a prophylactic therapy to prevent skin damage, e.g. pressure ulcers, and to reduce postoperative blistering.

Contra-indications

Arterial bleeds and heavily bleeding wounds.

Directions for use

Apply directly to the heel, pink side up and secure in place with tape, appropriate bandage or secondary dressing.

Frequency of change

For wound care, Advazorb® Heel can be left in place for up to 7 days but should be changed if exudate reaches 0.5cm from the edge of the dressing. Clinical observation is necessary to determine required frequency of change. For prophylactic therapy to prevent skin damage, inspect the skin regularly in line with local protocol and patient risk status. Carefully remove the dressing to allow assessment of skin and any bony prominence covered. Replace the dressing following assessment ensuring therapy is adjusted to take into account any changes. The dressing change interval may be several days or earlier if it is dislodged, rolled at the edges, wrinkled, creased or damaged, soiled, saturated or compromised in some other way.

Disposal

Dispose as clinical waste.

REF See Packaging	Catalogue #
	Single patient use
	Do not use if package is damaged
STERILE EO	Sterilized using ethylene oxide
	Do not re-sterilize
	Consult instructions for use
	Keep dry
	Keep away from sunlight
	Manufacturer
	Product is not made with natural rubber latex